


A photograph of a yellow concrete ledge with two metal handrails, set against a bright blue background. The handrails are made of dark metal and are mounted on the ledge. The concrete has some visible cracks and texture. The blue background is a solid, vibrant color.

STOCKS Training: Assay & Raw Data Management

Charles Girardot

Genome Biology Computational Support (GBCS)

From Specimen to Datasets: Complete Sample Lineage




Assays consume samples and generate raw datasets


- *Assay* is short for **raw data generating assay**
- Raw data is generated:
 - **from samples** prepared in a specific way e.g. a sequencing library
 - by **running a specific instrument** e.g. a sequencer or a microscope

- Raw data can be :

- One file
- Multiple files
- One folder
- Multiple folders
- A file (folder) is a *datafile*
- *Registering raw data == assay creation*



Assay types reflect data standard and user needs



- Assay types are specific to a particular experiment
- They belong to one of the technology:
 - Sequencing : Illumina Sequencing, Nanopore Sequencing
 - Light microscopy: LM or LM screen
 - [electron microscopy, proteomics, metabolomics, flow cytometry]
- Together with sample sub-types, they catch the FAIR aspects of a technology
 - allow easy data deposition
- Assay (and sample) sub-types also accommodate the user and institute needs
 - Custom properties available

The Illumina & Nanopore Sequencing Assay

STOCKS v1.26.5 (v1.21.10)

ASSAYS - NGSILLUMINAASSAY

HNHCMAX3_lane1

Name: HNHCMAX3_lane1
ID: 8f27b0da-b41c-4b08-a8d6-891c5d7853b2
State: REGISTERED
Multiplexed: yes Sample Number: 8
Workflow: —
Studies: [Cut&Tag_optimization_wing_discs](#) [nippedB_RNA](#)
Run Directory: [/g/furlong/STOCKS/Data/Assay/sequencing/2022/2022-06-02-HNHCMAX3](#)
Description: —

Run Parameters Sequencing specific

Run ID: [220531_NB501764_1566_AHNHCMAX3](#)
Instrument: [NextSeq500 Waffles](#)
Flowcell: HNHCMAX3 Flowcell Version: —
Lane: 1
Run Type: Paired-End
Run Mode: 75-8-8-75 Read Length: —
Chemistry: —

Sharing: —
Ownership and Lifecycle: —

Lineage Graph

Assays (1) | Datasets (8) | Samples (8)

8 samples: nlp0hrep2, nlp0hrep1, nlp0hrep2, nlp0hrep1, chipoldrep2, chipoldrep1, chipnewrep2, chipnewrep1

8 datasets: nlp0hrep1_HNHCMAX3_1, nlp0hrep2_HNHCMAX3_1, chipoldrep1_HNHCMAX3_1, chipoldrep2_HNHCMAX3_1, chipnewrep1_HNHCMAX3_1, chipnewrep2_HNHCMAX3_1

Output Datasets

Name	Description	Type	QC	Samples	Assay	Studies	Datafiles	Dataset Collections	Workflow	Archives
<input type="checkbox"/> Chipnewrep1_HNHCMAX3_1		paired-end fastq	—	Chipnewrep1	HNHCMAX3_lane1	Cut&Tag_optimization_wing_discs	view 2 items			
<input type="checkbox"/> Chipnewrep2_HNHCMAX3_1		paired-end fastq	—	Chipnewrep2	HNHCMAX3_lane1	Cut&Tag_optimization_wing_discs	view 2 items			
<input type="checkbox"/> Chipoldrep1_HNHCMAX3_1		paired-end fastq	—	Chipoldrep1	HNHCMAX3_lane1	Cut&Tag_optimization_wing_discs	view 2 items			

GBCS · 2018-present · [📧](#) [📧](#) [📧](#)

Assays have a multiplexing status (and sample number)

STOCKS v1.26.5 (v1.21.10)

ASSAYS - NGSILLUMINAASSAY

HNHCMAX3_lane1

Name: HNHCMAX3_lane1
ID: 8f27b0da-b41c-4b08-a8d6-891c5d7853b2
State: REGISTERED

Multiplexed ✓ yes **Sample Number** 8

Workflow: —
Studies: Cut&Tag_optimization_wing_discs, nippedB_RNA
Run Directory: /g/furlong/STOCKS/Data/Assay/sequencing/2022/2022-06-02-HNHCMAX3
Description: —

Run Parameters

Run ID: 220531_NB501764_1566_AHNHCMAX3
Instrument: NextSeq500 Waffles
Flowcell: HNHCMAX3 **Flowcell Version** —
Lane: 1
Run Type: Paired-End
Run Mode: 75-8-8-75 **Read Length** —
Chemistry: —

Sharing
Ownership and Lifecycle

Annotations 0
Notes 0

Lineage Graph

Assays (1) Datasets (8) Samples (8)

8 samples: nip0hrep2, nip0hrep1, nip0hrep2, nip0hrep1, chipoldrep2, chipoldrep1, chipnewrep2, chipnewrep1

8 datasets: nip0hrep1_HNHCMAX3_1, nip0hrep2_HNHCMAX3_1, chipoldrep1_HNHCMAX3_1, chipoldrep2_HNHCMAX3_1, chipnewrep1_HNHCMAX3_1, chipnewrep2_HNHCMAX3_1

Output Datasets

Name	Description	Type	QC	Samples	Assay	Studies	Datafiles	Dataset Collections	Workflow	Archives
Chipnewrep1_HNHCMAX3_1		paired-end fastq	—	Chipnewrep1	HNHCMAX3_lane1	Cut&Tag_optimization_wing_discs	view 2 items			
Chipnewrep2_HNHCMAX3_1		paired-end fastq	—	Chipnewrep2	HNHCMAX3_lane1	Cut&Tag_optimization_wing_discs	view 2 items			
Chipoldrep1_HNHCMAX3_1		paired-end fastq	—	Chipoldrep1	HNHCMAX3_lane1	Cut&Tag_optimization_wing_discs	view 2 items			

Also : Workflow can be linked ; Annotations and Notes as often

Easily navigate to consumed samples and produced datasets

The screenshot displays the STOCKS interface for an assay named "HNHCMF3_lane1". The interface is divided into several sections:

- Left Sidebar:** A navigation menu with categories like CONSUMABLES, EQUIPMENT, STORAGE EQUIPMENT, SPECIMEN, BIOMATERIALS, PROTOCOLS, LAB NOTEBOOK, ASSAYS, and DATASET MANAGEMENT. The "ASSAYS" section is currently active.
- Assay Details:** A central panel showing metadata for "HNHCMF3_lane1", including Name, ID, State, Multiplexed status, Sample Number (8), Workflow, Studies, Run Directory, and Description. Below this is a "Run Parameters" section with details like Run ID, Instrument (NextSeq500 Waffles), Flowcell, Lane, Run Type, Run Mode, and Chemistry.
- Annotations and Notes:** Two empty sections for adding annotations and notes to the assay.
- Lineage Graph:** A graph showing the relationship between the assay and its associated samples and datasets. The assay is at the center, connected to 8 samples (nlp0hrep1-2, chipnewrep1-2) and 8 datasets (nlp0hrep1_1-2, chipoldrep1-2, chipnewrep1-2).
- Output Datasets:** A table at the bottom listing the datasets produced by the assay. The table has columns for Name, Description, Type, QC, Samples, Assay, Studies, Datfiles, Dataset Collections, Workflow, and Archives.

Name	Description	Type	QC	Samples	Assay	Studies	Datfiles	Dataset Collections	Workflow	Archives
Chipnewrep1_HNHCMF3_1		paired-end fastq	—	Chipnewrep1	HNHCMF3_lane1	Cut&Tag_optimization_wing_discs	view 2 items			
Chipnewrep2_HNHCMF3_1		paired-end fastq	—	Chipnewrep2	HNHCMF3_lane1	Cut&Tag_optimization_wing_discs	view 2 items			

The assay datafiles are located in the *run directory*

STOCKS v1.26.5 (v1.21.10)

ASSAYS - NGSILLUMINAASSAY

HNHCMAFX3_lane1

Name: HNHCMAFX3_lane1
ID: 8f27b0da-b41c-4b08-a8d6-891c5d7853b2
State: REGISTERED
Multiplexed: yes Sample Number: 8
Workflow: —
Studies: [Cut&Tag_optimization_wing_discs](#), [nippedB_RNA](#)
Run Directory: /g/furlong/STOCKS/Data/Assay/sequencing/2022/2022-06-02-HNHCMAFX3
Description: —

Run Parameters

Run ID: [220531_NB501764_1566_AHHNHCMAFX3](#)
Instrument: [NextSeq500 Waffles](#)
Flowcell: HNHCMAFX3 Flowcell Version: —
Lane: 1
Run Type: Paired-End
Run Mode: 75-8-8-75 Read Length: —
Chemistry: —

Sharing

Ownership and Lifecycle

Annotations

Notes

Lineage Graph

Assays (1) | Datasets (8) | Samples (8)

8 samples: nip0hrep2, nip0hrep1, nip0hrep2, nip0hrep1, chipoldrep2, chipoldrep1, chipnewrep2, chipnewrep1

8 datasets: nip0hrep1_HNHCMAFX3_1, nip0hrep2_HNHCMAFX3_1, chipoldrep1_HNHCMAFX3_1, chipoldrep2_HNHCMAFX3_1, chipnewrep1_HNHCMAFX3_1, chipnewrep2_HNHCMAFX3_1

Output Datasets

Name	Description	Type	QC	Samples	Assay	Studies	Datafiles	Dataset Collections	Workflow	Archives
<input type="checkbox"/> Chipnewrep1_HNHCMAFX3_1		paired-end fastq	—	Chipnewrep1	HNHCMAFX3_lane1	Cut&Tag_optimization_wing_discs	view 2 items			
<input type="checkbox"/> Chipnewrep2_HNHCMAFX3_1		paired-end fastq	—	Chipnewrep2	HNHCMAFX3_lane1	Cut&Tag_optimization_wing_discs	view 2 items			
<input type="checkbox"/> Chipoldrep1_HNHCMAFX3_1		paired-end fastq	—	Chipoldrep1	HNHCMAFX3_lane1	Cut&Tag_optimization_wing_discs	view 2 items			

GBCS · 2018-present · [📧](#) [📧](#) [📧](#)

The Light Microscopy Assay Model

Pick a method, a microscope and an optional objective

Name	Golgi Morphology Imaging and Photoactivation Plate 1		
ID	dd0507b7-9f28-4ded-835b-7b5d05d517cb		
State	REGISTERED		
Imaging Method	confocal microscopy , fluorescence microscopy		
Instrument	ALMF-Zeiss LSM 900	—	
Objective	—		
Run Directory	/g/pepperkok/STOCKS/Data/Assay/light_microscopy/2022/2022_05_23_Zeiss_sasingh		
Description	The Neg9 cells were not imaged (directly photoactivated). This is why this imaging session does not contain Neg9 input sample nor images for it.		


Key values for the X, Y, Z and T dimensions

XYZCT Dimension Details					
Image X Size (px) <i>(-1 for varying)</i>	-1	Image Y Size (px) <i>(-1 for varying)</i>	-1	Pixel Size (microns)	-1
Z Stack	no	Z Planes	1	Z Step (microns)	-1
Time Lapse	no	Time Points (rounds)	1	Time Step (s)	0

Channels description (1+)

Channels				
#	Name	Target	Label	Probe
1	Autofocus	Nuclei (Nucleus)	DAPI (Fluorophore; Ex/Em max: 359/461 nm)	—
2	DAPI	Nuclei (Nucleus)	DAPI (Fluorophore; Ex/Em max: 359/461 nm)	—
3	Golgi	GalNac-GFP (Genetically encoded tag)	AcGFP1 (Fluorescent Protein; Ex/Em max: 475/505 nm)	—
4	Photoactivatable mCherry	H2B-PA-mCherry (Nucleus)	PAmCherry1(On) (Fluorescent Protein; Ex/Em max: 564/595 nm)	—

The Light Microscopy Assay Model – Channel description



Edit channel number 3

Channel Name
Golgi

Visualized Entity

Name* GalNac-GFP Type* Genetically encoded tag

Label

Name* AcGFP1 (Fluorescent Protein; Ex/Em max: 475/505 nm)

Object
Plasmid (Consumable)
Type to search Plasmid (Consumable)
dacB D3 His pET21b

Probe

Name Type

Object
Select type

Wavelength


Excitation Range (nm)
From To

Emission Range (nm)
From To

- Hundreds of Labels/Dyes imported from
 - Fluorescent proteins <https://www.fpbase.org/>
 - Fluorophores <https://www.semrock.com/>
- Also Designed with Flow Cytometry in mind
- Exact item used for the label & probe can be linked from the consumable collections

The Light Microscopy Assay Model – Microscopes


- All ALMF Microscopes available
- Groups can add their own microscopes



The screenshot shows the 'EQUIPMENT - MICROSCOPE' page in the STOCKS application. The left sidebar contains navigation options: 'CONSUMABLES', 'EQUIPMENT', 'Instrumentmodel', 'Store Items', 'Centrifuges', and 'Microscopes' (which is selected). The main content area displays a table of microscopes with columns for Name, Code, Origin/Vendor, Instrument Model, and Stock locations. The table lists several models, including Abbelight Safe360, Abberior STED/RESOLFT, Leica GSDIM, Leica SP8 STED 3x, Leica Stellaris 8, Leica TCS SP8 DLS, Luxendo Lightsheet InVi SPIM, and Luxendo Lightsheet MuVi SPIM. The 'Leica Microsystems' vendor name is highlighted in a tooltip.

Name	Code	Origin/Vendor	Instrument Model	Stock locations
ALMF-Abbelight Safe360	257	Abbelight	ALMF-Abbelight Safe360	13-509A-ALMF
ALMF-Abberior STED/RESOLFT	97	Abbelight	ALMF-Abberior STED/RESOLFT	13-516A-ALMF
ALMF-Leica GSDIM	98	Leica Microsystems	ALMF-Leica GSDIM	13-509D-ALMF
ALMF-Leica SP8 STED 3x	45	Leica Microsystems	ALMF-Leica SP8 STED 3x	13-516B-ALMF
ALMF-Leica Stellaris 8	313	Leica Microsystems	ALMF-Leica Stellaris 8	13-514D-ALMF
ALMF-Leica TCS SP8 DLS	51	Leica Microsystems	ALMF-Leica TCS SP8 DLS	13-509C-ALMF
ALMF-Luxendo Lightsheet InVi SPIM	157	Bruker	ALMF-Luxendo Lightsheet InVi SPIM	13-509E-ALMF
ALMF-Luxendo Lightsheet MuVi SPIM	89	Bruker	ALMF-Luxendo Lightsheet MuVi SPIM	13-509B-ALMF

- Code == PPMS ID
- Link to PPMS booking
- Responsible Person
- Room Location



The screenshot shows the details page for the 'ALMF-Abbelight Safe360' microscope. The page includes a header with 'EQUIPMENT - MICROSCOPE' and a search bar. Below the header, there are buttons for 'Clone', 'Copy URL', and 'Add to order'. The main content is divided into two columns. The left column, 'Item details', lists various attributes: ID (9b771de7-b875-4a8a-97b2-e5cc79fb46e0), Code (257), Origin/Vendor (Abbelight), Catalogue #, Description, Instrument Model (ALMF-Abbelight Safe360), Stock locations (1. 13-509A-ALMF), Status (AVAILABLE), Maintained (checked), Responsible Person (bneumann), Produced, Owner (bneumann), Owned by Group (Advanced Light Microscopy Core Facility), Modified (2022-04-28 11:55:00), Created (2022-04-28 11:55:00), and Permissions (Public item). The right column contains sections for 'Annotations' (1), 'Attachments' (0), 'Notes' (0), and 'Orders' (0). The 'Annotations' section shows a 'PPMS Link' with the URL 'https://ppms.embl.de/planning/?item=257'.

Item details	Value
ID	9b771de7-b875-4a8a-97b2-e5cc79fb46e0
Code	257
Origin/Vendor	Abbelight
Catalogue #	—
Description	—
Instrument Model	ALMF-Abbelight Safe360
Stock locations	1. 13-509A-ALMF
Status	AVAILABLE
Maintained	checked
Responsible Person	bneumann
Produced	—
Owner	bneumann
Owned by Group	Advanced Light Microscopy Core Facility
Modified	2022-04-28 11:55:00
Created	2022-04-28 11:55:00
Permissions	Public item

Multiplexed assays may contain datasets that belong to different users and projects

- An assay does not necessarily contain related datasets
- The assay, the produced and consumed samples can have different owners ; e.g.
 - The lab technicians performed the assay => they own the assay
 - Datasets (and samples) are owned by different researchers
- Related datasets, produced by one or more assays, are grouped into a **Study**

The screenshot displays two main sections: 'Ownership and Lifecycle' and 'Output Datasets'.

Ownership and Lifecycle

Owner	scholtal	Group	Furlong Group
Created	2014-04-04 00:00:00	By	admin
Modified	2022-09-23 15:57:47	By	admin
Deleted	✗ no		

Output Datasets

1 - 5 of 5

Owner	Name	Type	QC	Samples	Studies	Datafiles
monfort	mesoRNAseq_553-54_nuclei_lane6	paired-end fastq	FAILED	mesoRNAseq_553-54_nuclei	RNAseq (not cleaned-up)	view 2 items
monfort	mesoRNAseq_68h-Sa_lane6	paired-end fastq	PASSED	mesoRNAseq_68h-Sa	view 2 items	view 2 items
girardot	mesoRNAseq_46h-Sb_lane6	paired-end fastq	PASSED	mesoRNAseq_46h-Sb_2	view 2 items	view 2 items
girardot	mesoRNAseq_34h-Sb_lane6	paired-end fastq	PASSED	mesoRNAseq_34h-Sb_2	view 2 items	view 2 items
girardot	C3WEMACX:700999_6	paired-end fastq	—	view 4 items	view 2 items	view 2 items

The 'scholtal' owner name in the 'Ownership and Lifecycle' section and the 'Owner' column in the 'Output Datasets' table are highlighted with red boxes. A diagram in the top right shows a central node 'mesoRNAseq_C3WEMACX_6' with lines connecting to various data points.

Studies are at the basis of data publication to public repositories

The screenshot displays the 'STUDIES' interface. At the top, there's a search bar and navigation icons. Below that, the title of the study is 'ATAC-seq of wildtype or Dppa2 knockout ESC upon release of...'. The 'Item details' section includes fields for ID, Project, Design, PubMed ID, and Description. The 'Description' field contains a detailed paragraph about the experiment. Below this, there are fields for 'Nr. of datasets', 'Archives', 'Modified', 'Created', and 'Owned by Group'. A 'Permissions' table follows, showing access for groups and users. The 'Datasets' section at the bottom lists four datasets with columns for Name, Description, Type, QC, and Samples.

Item details	
ID	2735409c-3cc2-47a9-a7d7-799c7eb70695
Project	VC20_007 CnR in Dppa2KO
Design	biological replicate, genotype design
PubMed ID	
Description	Chromatin accessibility was assayed in <i>wildtype</i> or <i>Dppa2</i> knockout ESC after 26 days of release of the trigger imposed by epigenetic editing. Samples were collected in two clonal <i>knockout</i> and <i>wildtype</i> lines after sorting at FACS of cells which maintained a repressive Esg1-tdTomato (TOMneg) reporter expression after 26 days of DOX washout (release of the trigger).
Nr. of datasets	4
Archives	
Modified	2021-05-11 11:56:49
Modified by	girardot
Created	2021-05-05 14:34:19
Created by	carlini
Owner	carlini
Owned by Group	Hackett Group


Permissions	View	Edit	Delete
Groups			
Hackett Group	✓	✓	
Users			
carlini	✓	✓	✓
girardot	✓	✓	

Name	Description	Type	QC	Samples
KOB67dDWTOMneg_HG7WLBGXG_1		paired-end fastq	PASSED	KOB6_26dDW_TOI
KOC77dDWTOMneg_HG7WLBGXG_1		paired-end fastq	PASSED	KOC7_26dDW_TOI
WT126dDW_HG7WLBGXG_1		paired-end fastq	PASSED	WT1_26dDW
WT226dDW_HG7WLBGXG_1		paired-end fastq	PASSED	WT2_26dDW

- Related raw data must accumulate in a *study*
 - Produced at different times, in different runs
 - e.g. Chip-seq study
- Belongs to one project
- Datasets can belong to multiple studies
 - E.g. when re-used
- The Study is what gets submitted to data repository
- Your paper may have multiple studies of different data types
 - Different sequencing flavors e.g. rna-seq, chip-seq
 - Different technologies e.g. imaging, sequencing

Assay vs Study vs Project


- The **project** regroups different items i.e. samples & datasets (more in the future)
- The **study** regroups related datasets of the same experiment type, those are analyzed together
- The **assay** is the run of one instrument and regroup the produced datasets (==same run metadata)




A cross-technology sequencing/imaging “multi-omics” project

Visualizing sample - Batch3_HeLa_GalNaC

Assays (3) Datasets (201) Samples (6)



Data Ingestion Scenarios



Datasets are made of one or more Datafiles

STOCKS development v1.27.0 (v1.22.0) 147c39e

DATAFILES

+ New Search

1 - 10 of 138

Name	URI	Filetype	Datasets
<input type="checkbox"/> 10hr_Elav_rep1A_R1.fastq.gz	/g/furlong/STOCKS/Data/Assay/sequencing/2022/2022_08_15...	gz	10hr_Elav_rep1A
<input type="checkbox"/> 10hr_Elav_rep1A_R2.fastq.gz	/g/furlong/STOCKS/Data/Assay/sequencing/2022/2022_08_15...	gz	10hr_Elav_rep1A
<input type="checkbox"/> 10hr_Elav_rep1B_R1.fastq.gz	/g/furlong/STOCKS/Data/Assay/sequencing/2022/2022_08_15...	gz	10hr_Elav_rep1B
<input type="checkbox"/> 10hr_Elav_rep1B_R2.fastq.gz	/g/furlong/STOCKS/Data/Assay/sequencing/2022/2022_08_15...	gz	10hr_Elav_rep1B
<input type="checkbox"/> 10hr_Elav_rep2A_R1.fastq.gz	/g/furlong/STOCKS/Data/Assay/sequencing/2022/2022_08_15...	gz	
<input type="checkbox"/> 10hr_Elav_rep2A_R2.fastq.gz	/g/furlong/STOCKS/Data/Assay/sequencing/2022/2022_08_15...	gz	
<input type="checkbox"/> 10hr_Elav_rep2B_R1.fastq.gz	/g/furlong/STOCKS/Data/Assay/sequencing/2022/2022_08_15...	gz	
<input type="checkbox"/> 10hr_Elav_rep2B_R2.fastq.gz	/g/furlong/STOCKS/Data/Assay/sequencing/2022/2022_08_15...	gz	
<input type="checkbox"/> 10hr_Elav_rep3A_R1.fastq.gz	/g/furlong/STOCKS/Data/Assay/sequencing/2022/2022_08_15...	gz	
<input type="checkbox"/> 10hr_Elav_rep3A_R2.fastq.gz	/g/furlong/STOCKS/Data/Assay/sequencing/2022/2022_08_15...	gz	

DATASETS

3_TAG_RNASeq_RAL-303_10-12h_H0Y46BGXX_1

Item details

ID	6d5ff378-4f78-432f-9ece-d03aff4cf0c0
Description	—
Type	paired-end fastq
QC	PASSED
Samples	3_TAG_RNASeq_RAL-303_10-12h
Parent datasets	—
Assay	3_TAG_multi_RNA_Seq2_H0Y46BGXX_1
Studies	All TagSeq Datasets Drosophila melanogaster whole embryo RNA-Seq experiments during embryogenesis
Datafiles	read_1: 3_TAG_RNASeq_RAL-303_10-12h_LIB18703_RBA17928_1.txt.gz read_2: 3_TAG_RNASeq_RAL-303_10-12h_LIB18703_RBA17928_2.txt.gz

- One Paired-End FASTQ dataset is composed of 2 FASTQ *datafiles*
- Datasets can have as many datafiles as needed
- A dataset should be regarded as an atomic entity
 - i.e. considering a unique datafile would not make sense
- You always manipulate datasets in STOCKS

Flexible assay data import with the (raw) Dataset Loader

> Uploading **Raw data** *implies* to described the **Assay (e.g. Light Microscopy)** and the **parent sample(s)**

STOCKS v1.26.0 (v1.21.1)

ASSAY TYPE SELECTION

Loading Raw Data

Raw data is generated by an assay. Raw datasets are linked to (1) the assay and (2) the sample it derives from. Start by creating a **new assay**, either from scratch, or using a template or an assay as a template. Once the assay is created, you will be redirected to the dataset importer wizard

I want to add data to an existing assay ...

Import data ^

Import raw datasets

Import datasets

STORAGE EQUIPMENT

SPECIMEN

BIOMATERIALS

PROTOCOLS

LAB NOTEBOOK

ASSAYS

DATASET MANAGEMENT

Sequencing Assays

- Illumina Sequencing No Template Select assay Create new
- Nanopore Sequencing No Template Select assay Create new

Imaging Assays

- Light Microscopy Screen No Template Select assay Create new
- Light Microscopy No Template Select assay Create new

- Templates for frequent setups
- Re-use **your own assay as template**

=> The **Dataset Loader Wizard** drives the user through this process

DATASET LOADER

Select Data Build datasets Verify Assign samples

Continue

Data must be first organized into a “run folder”

- One folder == One assay
- Data (files, folders) must be in your user “dropbox” organized in a logical way
- Group data per assay as data for a unique assay can be loaded at once
- Individual files or complete folders can be selected as “data files”

The screenshot displays the 'DATASET LOADER' interface. At the top, there is a search bar and a menu. Below the search bar is a progress bar with four steps: 'Select Data' (active), 'Build datasets', 'Verify', and 'Assign samples'. A 'Continue' button is located to the right of the progress bar. The main area shows a file browser with a path: '/g/funcgen/base/stocks-development/dropbox/Computational_Support_Genome_Biology/girardot'. The file browser shows a list of files and folders. The 'PAIRED-END_DEMU_EXAMPLE/' folder is selected, and its contents are displayed in a table. All files in the table have a checked checkbox in the first column, indicating they are selected for upload.

Name	Created	Modified	Size
XYZ_PE-DEMU_01_1_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_01_2_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_02_1_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_02_2_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_03_1_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_03_2_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_04_1_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_04_2_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_05_1_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_05_2_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_06_1_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_06_2_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_07_1_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_07_2_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_08_1_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_08_2_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_09_1_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B
XYZ_PE-DEMU_09_2_sequence.txt.gz	March 29th, 2022	March 7th, 2022	0 B

Multi-file datasets can be assembled in the Dataset Builder

- Build multi-file datasets with the Dataset Builder
- Datasets can be organized in different dataset collections

DATASET LOADER Search... Menu

Select Data Build datasets Verify Assign samples Continue

Dataset Builder

 18 files selected

Collection #1 x + Create collection

Name: Collection #1

Number of datafile(s) per dataset: - 2 +

Raw data:

Assay type: NGSILLUMINAASSAY

Stage datasets to continue ↓

Dataset Name	Datafile 1	Datafile 2	
	Ngs_readtype: Read 1	Ngs_readtype: Read 2	<input checked="" type="checkbox"/> Stage all
	Filter Preset: Standard Genecore	Filter Preset: Standard Genecore	
	Filter: <input type="text" value="_1"/> RE Matching 9 items	Filter: <input type="text" value="_2"/> RE Matching 9 items	
1 XYZ_PE-DEMU_01_sequence.txt	XYZ_PE-DEMU_01_1_sequence.txt.gz	XYZ_PE-DEMU_01_2_sequence.txt.gz	<input checked="" type="checkbox"/> Stage
2 XYZ_PE-DEMU_02_sequence.txt	XYZ_PE-DEMU_02_1_sequence.txt.gz	XYZ_PE-DEMU_02_2_sequence.txt.gz	<input checked="" type="checkbox"/> Stage
3 XYZ_PE-DEMU_03_sequence.txt	XYZ_PE-DEMU_03_1_sequence.txt.gz	XYZ_PE-DEMU_03_2_sequence.txt.gz	<input checked="" type="checkbox"/> Stage
4 XYZ_PE-DEMU_04_sequence.txt	XYZ_PE-DEMU_04_1_sequence.txt.gz	XYZ_PE-DEMU_04_2_sequence.txt.gz	<input checked="" type="checkbox"/> Stage

Raw Datasets must be linked to Samples

- Link Datasets to new or existing samples, associate to study & owner

DATASET LOADER

Search... Menu

Select Data **Build datasets** **Verify** **Assign samples** Submit

[Load samples](#) [Load studies](#)

	Dataset	Collection	Sequencing Library	* Barcode	Files	* Study	* Owner
<input type="checkbox"/>	1 ✓ XYZ_PE-DEMU_01 New	Collection #1	XYZ_PE-DEMU_01 <small>Or load existing samples</small>	ACTGGC 1	...	201701_CaptureC_test x	girardot
<input type="checkbox"/>	2 ✗ XYZ_PE-DEMU_02 New	Collection #1	XYZ_PE-DEMU_02 <small>Or load existing samples</small>	<input type="text"/> <small>This field is mandatory.</small>	...	201701_CaptureC_test x	girardot
<input type="checkbox"/>	3 ✗ XYZ_PE-DEMU_03 New	Collection #1	XYZ_PE-DEMU_03 <small>Or load existing samples</small>	<input type="text"/> <small>This field is mandatory.</small>	...	201701_CaptureC_test x	girardot

- Once completed, the files and folders must be imported into STOCKS
 - Async process (can take quite some time)
 - A job is created
 - You get a notification once ready



Check the status of your async tasks

- Different tasks occur in an asynchronous fashion

The screenshot displays the 'STOCKS' application interface. On the left is a sidebar with navigation options: CONSUMABLES, EQUIPMENT, STORAGE EQUIPMENT, SPECIMEN, BIOMATERIALS, PROTOCOLS, and LAB NOTEBOOK. The main area is titled 'USER TASKS' and contains a table of tasks. A red box highlights the list view icon in the top right corner of the task list.




Task	Related Items	Category	Status	Triggered	Duration
Ingesting data files - 2022_08_15_11_46_47.35_girardot	Raw Fastq Novaseq_HTGTHDSX...	DATA MANAG	SUCCESS	a month ago	17m 12s
Syncing datasets to Galaxy	view 24 items	GALAXY DATA	SUCCESS	4 months ago	1m 16s
Freezing Experiment: RNA-seq Analysis for Khaloon (Ephrussi lab)	RNA-seq Analysis for Khaloon (...)	EXPERIMENT F	SUCCESS	8 months ago	11s
Syncing datasets to Galaxy	view 40 items	GALAXY DATA	SUCCESS	8 months ago	1m 14s
Freezing Experiment: Evaluation of Pareidolia Results	Evaluation of Pareidolia Results	EXPERIMENT F	SUCCESS	8 months ago	11s
Syncing datasets to Galaxy	view 40 items	GALAXY DATA	SUCCESS	9 months ago	1m 08s

Automated data ingestion using the STOCKS API



- No more data loss
- Information available at GeneCore is passed on
- Enforce to connect datasets with samples
- Enforce to collect minimum information
- User never manipulate the files
- Files remain accessible (read-only) for analysis
- Model to be extended to other Core Facilities

Files remain accessible but are managed by STOCKS



One STOCKS repo per group


- Read only to ensure traceability
- Each group pays for its storage
- Unix group sticky bit: group can read by default

Intuitive organization

User can import data from their “dropbox”

Dataset deletion rules differ between managed and unmanaged assays

- Automatically loaded assays are managed
 - They are automatically loaded from a trusted source
 - They must be complete (no post raw dataset addition) at load time
 - Uncomplete assays must not be loaded => get to your admin
 - They cannot be deleted by users
 - Admin can unload wrong assays
 - Raw datasets cannot be deleted unless a backup copy exist
- Manually loaded datasets and associated assays are unmanaged
 - User can delete datasets and assays as needed
 - You will soon be able to convert **unmanaged** assays/datasets to **managed** ones
 - Prevent mistakes to happen



EM Data Support coming soon

EM Sample Type

EM Sample preparation

Preparation method Freeze substitution
Filler

Grid Preparation Details (Room Temperature TEM)

Grid Model* Support film
Section thickness Immunogold (nm)

Grid Preparation Details (Cryo-TEM)

Grid Model*

Instrument Type Instrument
Gas mix Power (Wd)
Glow time (s)

Micropatterns?

Vitrification (Cryo-TEM)

Instrument or Model
Temperature (C) Humidity (%)
Blot paper type Cryogen
Sample volume (µl) Pulciacs used
Wait time (s) Blot time (s)
Blot force or Blot side

FIB Milling (Cryo-TEM)

Instrument Sputter coating (µd)
Milling angle (°) Sputter coating time (Sec)
Number of lamellae SIS coating time (Sec)
Rough milling currents (pA)
Fine milling current (pA) lamella width (µm)
Final lamella thickness (nm)

Array Tomography Details

Collection substrate Sample width (nm)
Amount of sections Section thickness (nm)
Coating material Coating thickness (nm)



Transmission EM Assay

Pre-Fill with template

Name

Imaging Method only offer TEM, HVEM, IVEM with TEM as default

TEM Type single-particle Tomography

Assay Temperature Cryo Room Temperature

Electron Microscope

Description

Workflow/Protocol

Dataset Acquisition details

Pixel Size (Å)	<input type="text"/>	Nominal Magnification (X)	<input type="text"/>
C2 Aperture (microns)	<input type="text"/>	Objective Aperture (microns)	<input type="text"/>
Exposure Time (s)	<input type="text"/>	Electron Dose (e/A2)	<input type="text"/>
Image Format	<input type="text"/>	Movie Fraction	<input type="text"/>
Beam Size (µm)	<input type="text"/>	Spot Size	<input type="text"/>
Beam Mode	<input type="text"/>	Energy Filter	<input type="checkbox"/>
Slit width (eV)	<input type="text"/>	Montage	<input type="checkbox"/>

Tomography details:

Tilting Method

Series Size Angle step (degrees)

Tilt Range (degree): From To

or

Volume EM (Assay)

Pre-Fill with template

Name

Imaging Method FIB-SEM, SBFSEM and Array Tomography (no default)

Is cryo FIB-SEM? Only visible if FIB-SEM selected

Electron Microscope

Description

Workflow/Protocol

Dataset Acquisition details

Detector type	<input type="text"/>	Pixel Size (x/y)	<input type="text"/> <input type="text"/>
High-tension (kV)	<input type="text"/>	Primary beam current (pA)	<input type="text"/>
Image/Tile X Size (px)	<input type="text"/>	Image/Tile y Size (px)	<input type="text"/>
Tile overlap (SBEM)	<input type="text"/>	Amount of images	<input type="text"/>
Dwell time (µs)	<input type="text"/>	Line average	<input type="text"/>
Slice thickness/Z-spacing (nm)	<input type="text"/>		
FIB-SEM Milling current (pA)	<input type="text"/>		

=> Flow Cytometry and other omics assays to follow !

Thank you for your attention



Jelle Scholtalbers



Matthias Monfort



Nayeem Reza

Questions ?