


A photograph of a yellow concrete ledge with two metal handrails, set against a bright blue background. The handrails are made of dark metal and are mounted on the ledge. The concrete has some visible cracks and texture. The blue background is a solid, vibrant color.

STOCKS Training: Derived Datasets and Analysis

Charles Girardot

Genome Biology Computational Support (GBCS)

STOCKS aims at storing the analyzed datasets



Raw and derived datasets are listed under the same menu

The screenshot displays the 'DATASETS' management interface. On the left, a sidebar contains navigation options: PROTOCOLS, LAB NOTEBOOK, ASSAYS, DATASET MANAGEMENT (selected), Study, Dataset Collections, Dataset (highlighted), DataFile, and Archive. The main panel shows a table of datasets with the following columns: Modified, Name, Assay, Studies, Owner, QC, and Samples. The table contains 10 rows of data, including entries like 'plate2well1', 'metadata.txt', and 'XYZ_PE-DEMU_01_1_sequence.txt'. Search filters are visible above the table, and the interface includes a top navigation bar with a user profile '@admin'.

Modified	Name	Assay	Studies	Owner	QC	Samples
2022-09-23 17:22:19	plate2well1	Awesome Sequencing	MyStudy	scholtal	PASSED	plate2well1
2022-09-23 17:22:19	metadata.txt	Awesome Sequencing	MyStudy	scholtal	PASSED	metadata.txt
2022-09-23 17:22:19	XYZ_PE-DEMU_01_1_sequence.txt	Awesome Sequencing	MyStudy	scholtal	PASSED	XYZ_PE-DEMU_01_1_sequence.txt
2022-09-23 17:22:19	plate1_well1-VF1.czi		MyStudy	scholtal	PASSED	
2022-09-23 17:22:19	movies	Awesome Sequencing	MyStudy	scholtal	PASSED	movies
2022-09-23 17:22:19	metadata.txt		MyStudy	scholtal	PASSED	
2022-09-23 17:22:19	XYZ_PE-DEMU_01_1_sequence.txt	Awesome Sequencing	MyStudy	scholtal	PASSED	XYZ_PE-DEMU_01_1_sequence.txt
2022-09-23 17:22:19	movies		MyStudy	scholtal	PASSED	
2022-09-23 17:22:19	plate1_well1-VF1.czi		MyStudy	scholtal	PASSED	
2022-09-23 17:22:19	MyStudy	scholtal	PASSED	...

Raw data processing made easy with Galaxy Link

ATAC-seq of wildtype or Dppa2 knockout ESC upon release of...

Item details

ID: 2735409c-3cc2-47a9-a7d7-799c7eb70695

Project: VC20_007 CnR in Dppa2KO

Design: biological replicate, genotype design

PubMed ID:

Description: Chromatin accessibility was assayed in *wildtype* or *Dppa2* knockout ESC after 26 days of release of the trigger imposed by epigenetic editing. Samples were collected in two clonal *knockout* and *wildtype* lines after sorting at FACS of cells which maintained a repressive Esg1-tdTomato (TOMneg) reporter expression after 26 days of DOX washout (release of the trigger).

Nr. of datasets: 4

Archives:

Modified: 2021-05-11 11:56:49

Modified by: girardot

Created: 2021-05-05 14:34:19

Created by: carlini

Owner: carlini

Owned by Group: Hackett Group

Permissions:

	View	Edit	Delete
Groups			
Hackett Group	✓	✓	
Users			
carlini	✓	✓	✓
girardot	✓	✓	

Datasets

Name	Description	Type	QC	Samples
KOB67dWTOMneg_HG7WLBGXG_1		paired-end fastq	PASSED	KOB6_26dW_TOI
KOC77dWTOMneg_HG7WLBGXG_1		paired-end fastq	PASSED	KOC7_26dW_TOI
WT126dDW_HG7WLBGXG_1		paired-end fastq	PASSED	WT1_26dDW
WT226dDW_HG7WLBGXG_1		paired-end fastq	PASSED	WT2_26dDW

No data duplication →


<https://galaxy.embl.de>

Galaxy

Libraries / Hackett Lab / VC20_007 ATAC-seq Dppa2KO

Name	Description	Type	Size	Updated	State
HG7WLBGXG_VC20_006-VC20_007_20x002698-1-1_Carlini_lane1KOB67dWTOMneg_1_sequence.txt.gz	uploaded fastqsanger.gz file	fastqsanger.gz	845.3 MB	a year ago	
HG7WLBGXG_VC20_006-VC20_007_20x002698-1-1_Carlini_lane1KOB67dWTOMneg_2_sequence.txt.gz	uploaded fastqsanger.gz file	fastqsanger.gz	842.1 MB	a year ago	
HG7WLBGXG_VC20_006-VC20_007_20x002698-1-1_Carlini_lane1KOB67dWTOMneg_1_sequence.txt.gz	uploaded fastqsanger.gz file	fastqsanger.gz	241.6 MB	a year ago	
HG7WLBGXG_VC20_006-VC20_007_20x002698-1-1_Carlini_lane1KOB67dWTOMneg_2_sequence.txt.gz	uploaded fastqsanger.gz file	fastqsanger.gz	242.5 MB	a year ago	
HG7WLBGXG_VC20_006-VC20_007_20x002698-1-1_Carlini_lane1KOB67dWTOMneg_1_sequence.txt.gz	uploaded fastqsanger.gz file	fastqsanger.gz	497.6 MB	a year ago	
HG7WLBGXG_VC20_006-VC20_007_20x002698-1-1_Carlini_lane1KOB67dWTOMneg_2_sequence.txt.gz	uploaded fastqsanger.gz file	fastqsanger.gz	498.1 MB	a year ago	
HG7WLBGXG_VC20_006-VC20_007_20x002698-1-1_Carlini_lane1KOC77dWTOMneg_1_sequence.txt.gz	uploaded fastqsanger.gz file	fastqsanger.gz	419.4 MB	a year ago	
HG7WLBGXG_VC20_006-VC20_007_20x002698-1-1_Carlini_lane1KOC77dWTOMneg_2_sequence.txt.gz	uploaded fastqsanger.gz file	fastqsanger.gz	421.1 MB	a year ago	
HG7WLBGXG_VC20_006-VC20_007_20x002698-1-1_Carlini_lane1KOC77dWTOMneg_1_sequence.txt.gz	uploaded fastqsanger.gz file	fastqsanger.gz	1.1 GB	a year ago	
HG7WLBGXG_VC20_006-VC20_007_20x002698-1-1_Carlini_lane1KOC77dWTOMneg_2_sequence.txt.gz	uploaded fastqsanger.gz file	fastqsanger.gz	1.1 GB	a year ago	

Start ATAC-seq Standard WF



All jobs executed on ITS HPC

Analyzed Data

Sending derived datasets back to STOCKS


Under active dev

The screenshot displays the Galaxy interface. On the left, a 'History' panel shows a list of datasets, with the entry '36: Send to STOCKS on data 32 and data 35' highlighted in orange. A red arrow points from this entry to the main workflow view. The main view shows the details for the workflow step 'Paste on data 46 and data 45'. The 'Item details' section includes fields for ID, Description, Type, QC, Samples, Assay, Studies, and Datafiles. The 'Dataset Collections' section is highlighted with a red box, showing 'Workflow' as 'Galaxy Workflow: Concat'. The 'Annotations' section shows the Galaxy History Name and ID. The 'Lineage' section, also highlighted with a red box, shows a graph of the workflow steps: 'DSM 17678' (specimen) leads to 'D1' (sample), which leads to 'Concatenate datasets on data 2 and data 1', 'Remove beginning on data 44', and 'Convert on data 44'. The 'Paste on data 46 and data 45' step is the final output of this lineage.


- Derived datasets are **automatically** linked to their parent
- The **complete data structure** is reproduced in STOCKS
- The executed **WF is saved as a protocol**
- **Links** to Galaxy are maintained

Processed datasets can also be manually imported

From the main Menu



From the Assay Page




- Loading derived datasets uses the same dataset loader, but simplified
 - No Assay type selection nor Assay creation
 - Datasets can optionally linked to individual samples
- Linking derived datasets to samples is highly recommended (when relevant)



Import multiple Dataset Collection at once

- STOCKS also offers a **DatasetCollection** to (manually) group datasets into meaningful sets
 - e.g. “bam files”, “segmentation maps”
- They are created when importing datasets



The Dataset Builder interface shows a collection named 'Aligned reads' with 2 collections. The 'Advanced collection controls' panel is open, showing the 'Dataset name extractor' field with the regex pattern `^([A-z]+)_*([0-9]+)*$`. The main table lists datasets with their names and datafile paths.

Dataset Name	Datafile 1	Filter	Stage
1 XYZ_01	XYZ_01_nodups.bam	e.g. _1.fastq.gz	Stage
2 XYZ_02	XYZ_02_nodups.bam		Stage
3 XYZ_03	XYZ_03_nodups.bam		Stage
4 XYZ_04	XYZ_04_nodups.bam		Stage

Linking processed datasets to parent datasets

DATASETS

movies

Item details

ID * 197ac00b-3366-4907-b2ec-1749d26e1050

Description

Type * generic

QC PASSED


Samples —

Parent datasets


Type to search Dataset

		1. Fragmented Plate3
		2. Fragmented plate2
		3. Fragmented plate1


Assay —



STOCKS to support complex Multi-omics scenarios



- Aggregate on common samples
- Dependent Datasets can be linked



Datafiles may have multiple copies (different locations)

- A paired-end datasets (== 2 datafiles) with two copies (hence 4 datafile copies)

The screenshot shows a web interface for managing datasets. The top header is green and contains the word "DATASETS" and a user profile "@admin". Below the header, the dataset name "3_TAG_RNASeq_RAL-303_10-12h_H0Y46BGXX_1" is displayed with "Edit" and "Copy URL" buttons. A "Groups" section shows "Furlong Group" with a checkmark. A "Users" section shows "cannavo" with checkmarks in three columns. Below this is a "Datafile copies" section with a search bar and filters for "All", "Group", and "Personal". A table lists 4 datafile copies with columns for URI, Assays, Short name, Peek, and Datafile Version.

URI	Assays	Short name	Peek	Datafile Version
/filerarchive/gbcs/embase_archive_51-girardot-Variou...	3_TAG_multi_RNA_Seq2_H0Y46BGXX_1	3_TAG_RNASeq_RAL-303_	—	1: 3_TAG_RNASeq_RAL-303_10-12h_LIB18703_RBA17928_1.txt...
/filerarchive/gbcs/embase_archive_51-girardot-Variou...	3_TAG_multi_RNA_Seq2_H0Y46BGXX_1	3_TAG_RNASeq_RAL-303_	—	1: 3_TAG_RNASeq_RAL-303_10-12h_LIB18703_RBA17928_2.txt...
/g/furlong/STOCKS/Data/Assay/sequencing/2015/2015-01-08-...	3_TAG_multi_RNA_Seq2_H0Y46BGXX_1	3_TAG_RNASeq_RAL-303_	—	1: 3_TAG_RNASeq_RAL-303_10-12h_LIB18703_RBA17928_1.txt...
/g/furlong/STOCKS/Data/Assay/sequencing/2015/2015-01-08-...	3_TAG_multi_RNA_Seq2_H0Y46BGXX_1	3_TAG_RNASeq_RAL-303_	—	1: 3_TAG_RNASeq_RAL-303_10-12h_LIB18703_RBA17928_2.txt...

- Tracking copies is important:
 - to save space by deleting useless copies
 - All copies of human data must be deleted when project ends
 - STOCKS has a concept of **primary copy** (cannot be deleted but in special situations)

Data can be archived to free space

The screenshot shows the STOCKS ARCHIVES interface. The top navigation bar is green with the text 'ARCHIVES' and a user profile '@admin'. Below the navigation bar is a search bar and a dropdown menu for 'All', 'Group', and 'Personal'. A pagination bar shows '1 - 10 of 29' items. The main content area is a table with columns: Name, URI, Projects, Studies, and Nr. of datasets. The table lists several archives, including '3Tfill_Wilkening_DKFZ', '3Tfil_YJR141W_woth_charlie_Boone_lab', '768_genotypes_experiment_archive', 'All_Avner_Data_01_June_2020', 'All_CAGE_PRO-cap_and_lncRNA_from_schor_and_males', 'All NGS Data till July 2020 - Ephrussi Group', 'Archive_2014-07_2016-01', 'Archive_5TagSeq_fromGrave', 'CB_index_failed', and 'Clouder Archiving Request March2020'. A sidebar on the left contains navigation options: 'Import datasets', 'LAB NOTEBOOK', 'ASSAYS', 'DATASET MANAGEMENT', 'Study', 'Dataset Collections', 'Dataset', 'DataFile', and 'Archive' (highlighted).

Name	URI	Projects	Studies	Nr. of datasets
3Tfill_Wilkening_DKFZ	/filerarchive/gbcs/			15
3Tfil_YJR141W_woth_charlie_Boone_lab	/filerarchive/gbcs/			8
768_genotypes_experiment_archive	/filerarchive/gbcs/			18
All_Avner_Data_01_June_2020	/filerarchive/gbcs/			492
All_CAGE_PRO-cap_and_lncRNA_from_schor_and_males	/filerarchive/gbcs/			468
All NGS Data till July 2020 - Ephrussi Group	/filerar	NGS data till July 2020 - Ephrussi Group	view 56 items	129
Archive_2014-07_2016-01	/filerarchive/gbcs/			
Archive_5TagSeq_fromGrave	/filerarchive/gbcs/			
CB_index_failed	/filerarchive/gbcs/			
Clouder Archiving Request March2020	/filerarchive/gbcs/			

The screenshot shows the details for the 'Archive_2014-07_2016-01' archive. The top navigation bar is green with the text 'ARCHIVES' and a user profile '@admin'. Below the navigation bar is a search bar and a dropdown menu for 'All', 'Group', and 'Personal'. A pagination bar shows '1 - 10 of 114' items. The main content area is divided into two sections: 'Datasets' and 'Datafile copies'. The 'Datasets' section has a search bar and a dropdown menu for 'All', 'Group', and 'Personal'. It shows a table with columns: Name, Description, Type, QC, Samples, Parent datasets, and Assay. The 'Datafile copies' section has a search bar and a dropdown menu for 'All', 'Group', and 'Personal'. It shows a table with columns: URI, Assays, Short name, Peek, and Datafile Version.

Name	Description	Type	QC	Samples	Parent datasets	Assay
3T_cytRNA-Aniso-RII_lane1		paired-end fastq	NA	3T_cytRNA-Aniso-RII		3T_C2R8HA
3T_cytRNA-Aniso-R-I_lane1		paired-end fastq	NA	3T_cytRNA-Aniso-R-I		3T_C2R8HA
3T_cytRNA-DMSO-RIII_lane1		paired-end fastq	NA	3T_cytRNA-DMSO-RIII		3T_C2R8HA
3T_cytRNA-DMSO-R-II_lane1		paired-end fastq	NA	3T_cytRNA-DMSO-R-II		3T_C2R8HA
3T_cytRNA-DMSO-R-I_lane1		paired-end fastq	NA	3T_cytRNA-DMSO-R-I		3T_C2R8HA
AM-GFP_lane8		paired-end fastq	NA	AM-GFP		C4VKHAC0
AM-parental-1_lane8		paired-end fastq	NA	AM-parental-1		C4VKHAC0
AM-parental-2_lane8		paired-end fastq	NA	AM-parental-2		C4VKHAC0
AM-parental-3_lane8		paired-end fastq	NA	AM-parental-3		C4VKHAC0

URI	Assays	Short name	Peek	Datafile Version
/filerarchive/gbcs/embase_archive_46-biohentze-Archive_201...	MOV10_C258FACXX_8	C258FACXX_JCLIPMOV10_	—	1: C258FACXX_JCLIPMOV10_13s0C
/filerarchive/gbcs/embase_archive_46-biohentze-Archive_201...	hnRNPC_C281TACXX_8	C281TACXX_JCLIPhnRNPC	—	1: C281TACXX_JCLIPhnRNPC_13s0C

- Old data is packed in an **archive and copied to Tape**
- All **metadata** remains live and searchable
- Ultimately handled by ITS
- Taping cost apply (see ITS)

Smart Dataset Archiving

STOCKS development v1.25.2 (v1.20.4) dfbf58

STUDIES

1 - 29 of 29 (5 selected)

<input type="checkbox"/>	Nr. of datasets	Name	Array	Project
<input type="checkbox"/>	0	my new study on the fly	—	Tra
<input checked="" type="checkbox"/>	10	Yad Hi-C 6-8h (2017)	—	Top
<input checked="" type="checkbox"/>	10	Yad Hi-C 2-4h (2017)	—	em
<input checked="" type="checkbox"/>	10	Yad HiC 16-18h (2017)	—	Top
<input checked="" type="checkbox"/>	22	BITSHeart_ChIPseq_all	—	em
<input checked="" type="checkbox"/>	12	BiTSHeart FAIRE	—	em
<input type="checkbox"/>	14	BiTS_4-6h	—	em
<input type="checkbox"/>	0	Toll10B_MZT	—	em
<input type="checkbox"/>	0	BiTS_6-8h	—	em
<input type="checkbox"/>	0	BiTS_2-4h	—	em
<input type="checkbox"/>	0	Drosophila melanogaster strain sequencing	—	em
<input type="checkbox"/>	0	MedakaChIP	—	em
<input type="checkbox"/>	0	BiTS_K9me2-3	—	em
<input type="checkbox"/>	0	Nau_ChIP-seq	—	em
<input type="checkbox"/>	0	MZT_GROSeq	—	em

Archive selected items

Batch Archive Items

5 Studies

62 Datasets (92 datafiles selected, 60 archivable) 32 datafiles cannot be archived and will be ignored

Dataset	Datafile	Size	Status	Archivable
BITSHeart_input1_lane2	BITSHeart_input1_[Heart1]_K4me3_10-13h_LIB16293_RBA15954_1.txt.gz	0 B		✗
BITSHeart_FAIRE4_lane2	BITSHeart_FAIRE4_10-13h_LIB16294_RBA15955_1.txt.gz	0 B		✗
BITSHeart_FAIRE6_lane2	BITSHeart_FAIRE6_13-16h_LIB16295_RBA15956_1.txt.gz	0 B		✗
BITSHeart_input2_lane2	BITSHeart_input2_10-13h_LIB16296_RBA15957_1.txt.gz	0 B		✗
BITSHeart_input4_lane2	BITSHeart_input4_[Heart14]_K9me3_13-16h_LIB16297_RBA15958_1.txt.gz	0 B		✗
BITSHeart_Heart1_lane2	BITSHeart_Heart1_K4me3_10-13h_LIB16298_RBA15959_1.txt.gz	0 B		✗
Total Size		675 GB		

Archive Information (* indicates a required field)

Name * Billing Details *

Description

In this archive, we have 60 datafiles of size - 675 GB from the following

Studies:

- Yad Hi-C 6-8h (2017) - gbcs-dev.embl.de/stocks/master/304e679f-d4c8-46a5-b325-cf3a63a21846

Ignore datafiles which are already backed-up/scheduled for a backup

1. Tape space will be billed to your group
2. Archiving is an asynchronous operation that can take up to a few days to complete
3. Once the archive is finished, deletion options will be available

I understand and confirm the above information is being cross-checked! *

Archive

Close

- Send Datasets to tape by **selecting Studies or Projects**
- Datasets have a **new copy** on tape
- Datasets can be deleted from disks **after archiving**
 - This removes the main copy of the datasets
 - The **tape copy is immutable** and cannot be deleted; kept for 10 years
- Already archived can be skipped

Datafiles are always listed, even when offline

The screenshot displays the STOCKS DataFiles interface. The left sidebar contains navigation options: ASSAYS, DATASET MANAGEMENT (with sub-items: Study, Dataset Collections, Dataset, DataFile, Archive), and GBCS (with sub-items: 2018-present, and icons for a document, chat, and email). The main content area is titled 'DATAFILES' and shows a list of datafiles. The list has columns for ID, Name, URI, Filetype, Datasets, and Deleted. A search bar at the top left of the list contains the text 'mnase'. The list shows 16 items, with the first four items highlighted in a red box. These four items have the status 'No file copy online' in the URI column. The interface also includes a '+ New' button, a search bar, and pagination controls showing '1 - 10 of 16'.

ID	Name	URI	Filetype	Datasets	Deleted
<input type="checkbox"/>	0d944013	000000000-A2839_MNaseEDTA_13s000141-1-1_Chabbert_lan... No file copy online	fastq	000000000-A2839:MiSeqSe01_1	<input checked="" type="checkbox"/> no
<input type="checkbox"/>	26a6228c	000000000-A2839_MNaseEDTA_13s000141-1-1_Chabbert_lan... No file copy online	fastq	000000000-A2839:MiSeqSe01_1	<input checked="" type="checkbox"/> no
<input type="checkbox"/>	2ea72270	A2EKD_MnaseLIB121912_12s008126-1-1_Adjalley_lane112s00... No file copy online	fastq	A2EKD:MiSeqSe01_1	<input checked="" type="checkbox"/> no
<input type="checkbox"/>	99feade5	A2EKD_MnaseLIB121912_12s008126-1-1_Adjalley_lane112s00... No file copy online	fastq	A2EKD:MiSeqSe01_1	<input checked="" type="checkbox"/> no
<input type="checkbox"/>	9bc82948	C6NAKACXX_TIGR_MNase_15s009387-1-1_Mikhaylichenko_Jan... /g/furlong/STOCKS/Data/Assay/sequencing/2015/2015-03-13-...	fastq	C6NAKACXX:700169_6	<input checked="" type="checkbox"/> no
<input type="checkbox"/>	5e9e54a4	C8Y21ACXX_RAI1-PHF14-K27ac_16s004082-1-1_Serandour_la... /g/noh/STOCKS/Data/Assay/sequencing/2016/2016-06-28-C8Y...	fastq	ChIPMnaseH3K27acD0_C8Y21ACXX_lane3	<input checked="" type="checkbox"/> no
<input type="checkbox"/>	1a57f813	C8Y21ACXX_RAI1-PHF14-K27ac_16s004082-1-1_Serandour_la... /g/noh/STOCKS/Data/Assay/sequencing/2016/2016-06-28-C8Y...	fastq	ChIPMnaseH3K27acD8_C8Y21ACXX_lane3	<input checked="" type="checkbox"/> no
<input type="checkbox"/>	76366445	HGLLLBGXX_gcn5_MNase_16s000049-1-1_Schiklenk_lane1134... /g/haering/ngsdata_lib/2016-01-19-HGLLLBGXX/lane1/LIB209...	fastq	gcn5MN_LY1349N1_HGLLLBGXX_lane1	<input checked="" type="checkbox"/> no
<input type="checkbox"/>	3c1323dd	HGLLLBGXX_gcn5_MNase_16s000049-1-1_Schiklenk_lane1134... /g/haering/ngsdata_lib/2016-01-19-HGLLLBGXX/lane1/LIB209...	fastq	gcn5MN_LY1349N2_HGLLLBGXX_lane1	<input checked="" type="checkbox"/> no
<input type="checkbox"/>	5b400168	HGLLLBGXX_gcn5_MNase_16s000049-1-1_Schiklenk_lane1134... /g/haering/ngsdata_lib/2016-01-19-HGLLLBGXX/lane1/LIB209...	fastq	gcn5MN_LY1349N3_HGLLLBGXX_lane1	<input checked="" type="checkbox"/> no

Preparing data for public repository submission

- Data submission == Study submission
- Make sure to have a study containing all the datasets to submit
 - Give a good name and description (will be public)
 - Set the design terms; which also hint at necessary sample annotations
- Annotate your samples according to the study design
- Make sure all samples' relevant fields are filled
- Link protocols to your samples (all the way up to specimen)
 - Make sure each protocol has a **summary** (no format, few sentences)
 - Growth protocol is always mandatory
 - Extraction protocol is always mandatory for omics
 - NGS Library preparation is always mandatory for sequencing
- One ready get in touch with GBCS
- Get in touch a few weeks before you want to submit your paper
 - Get in touch at project start if you use human samples

Data Submission as done today

Usual submission steps (after initial discussion):

- Annotate samples
- Curate samples (e.g. merging redundant samples)
- Linking protocols
- *Visualization of the complete Study Graph*



STOCKS UI

- Export MAGETAB
 - \$ stocks arrayexpress-export export.xlsx <studyid>
- Transfer Files
 - Create local links to files
 - FTP or Aspera
- Email the MAGE-TAB to annotate



STOCKS CLI (uses REST API)

<https://git.embl.de/grp-gbcs/stocks-python-client>


```
(/g/funcgen/galaxy/conda-envs/stocks-client) galaxy@spinoza:/g/funcgen/embase_submissions/kaesper_april2021/rai1/fastq$ stocks list datafilecopies 6ae41773-1245-4e12-95a3-1ae9379a26a5 --query_param "page_size=max" | jq -r '.results[]|[.uri] | @tsv' | awk '{system("ln -s \"$1\"")}'
Filter by UUID: 6ae41773-1245-4e12-95a3-1ae9379a26a5 (study)
(/g/funcgen/galaxy/conda-envs/stocks-client) galaxy@spinoza:/g/funcgen/embase_submissions/kaesper_april2021/rai1/fastq$ ls -l
total 48
lrwxrwxrwx 1 galaxy gbcs 156 Apr 28 15:24 CBK33ACXX_RNAD12-KORAI1ex4_16s001173-1-1_Serandour_lane6RNAES129D12c11_sequence.txt.gz -> /g/noh/NGS_Library/2016-03-10-CBK33ACXX/lane6/LIB21152_RBA20535/fastq/CBK33ACXX_RNAD12-KORAI1ex4_16s001173-1-1_Serandour_lane6RNAES129D12c11_sequence.txt.gz
lrwxrwxrwx 1 galaxy gbcs 156 Apr 28 15:24 CBK33ACXX_RNAD12-KORAI1ex4_16s001173-1-1_Serandour_lane6RNAES129D12c12_sequence.txt.gz -> /g/noh/NGS_Library/2016-03-10-CBK33ACXX/lane6/LIB21153_RBA20536/fastq/CBK33ACXX_RNAD12-KORAI1ex4_16s001173-1-1_Serandour_lane6RNAES129D12c12_sequence.txt.gz
lrwxrwxrwx 1 galaxy gbcs 156 Apr 28 15:24 CBK33ACXX_RNAD12-KORAI1ex4_16s001173-1-1_Serandour_lane6RNAES129D12c14_sequence.txt.gz -> /g/noh/NGS_Library/2016-03-10-CBK33ACXX/lane6/LIB21151_RBA20534/fastq/CBK33ACXX_RNAD12-KORAI1ex4_16s001173-1-1_Serandour_lane6RNAES129D12c14_sequence.txt.gz
lrwxrwxrwx 1 galaxy gbcs 156 Apr 28 15:24 CBK33ACXX_RNAD12-KORAI1ex4_16s001173-1-1_Serandour_lane6RNAES129D12c15_sequence.txt.gz -> /g/noh/NGS_Library/2016-03-10-CBK33ACXX/lane6/LIB21149_RBA20532/fastq/CBK33ACXX_RNAD12-KORAI1ex4_16s001173-1-1_Serandour_lane6RNAES129D12c15_sequence.txt.gz
lrwxrwxrwx 1 galaxy gbcs 156 Apr 28 15:24 CBK33ACXX_RNAD12-KORAI1ex4_16s001173-1-1_Serandour_lane6RNAES129D12c16_sequence.txt.gz -> /g/noh/NGS_Library/2016-03-10-CBK33ACXX/lane6/LIB21150_RBA20533/fastq/CBK33ACXX_RNAD12-KORAI1ex4_16s001173-1-1_Serandour_lane6RNAES129D12c16_sequence.txt.gz
lrwxrwxrwx 1 galaxy gbcs 156 Apr 28 15:24 CBK33ACXX_RNAD12-KORAI1ex4_16s001173-1-1_Serandour_lane6RNAES129D12c18_sequence.txt.gz -> /g/noh/NGS_Library/2016-03-10-CBK33ACXX/lane6/LIB21148_RBA20531/fastq/CBK33ACXX_RNAD12-KORAI1ex4_16s001173-1-1_Serandour_lane6RNAES129D12c18_sequence.txt.gz
lrwxrwxrwx 1 galaxy gbcs 154 Apr 28 15:24 CBK33ACXX_RNAD8-KORAI1ex4_16s001174-1-1_Serandour_lane7RNAES129D8c11_sequence.txt.gz -> /g/noh/NGS_Library/2016-03-10-CBK33ACXX/lane7/LIB21146_RBA20529/fastq/CBK33ACXX_RNAD8-KORAI1ex4_16s001174-1-1_Serandour_lane7RNAES129D8c11_sequence.txt.gz
lrwxrwxrwx 1 galaxy gbcs 154 Apr 28 15:24 CBK33ACXX_RNAD8-KORAI1ex4_16s001174-1-1_Serandour_lane7RNAES129D8c12_sequence.txt.gz -> /g/noh/NGS_Library/2016-03-10-CBK33ACXX/lane7/LIB21144_RBA20527/fastq/CBK33ACXX_RNAD8-KORAI1ex4_16s001174-1-1_Serandour_lane7RNAES129D8c12_sequence.txt.gz
lrwxrwxrwx 1 galaxy gbcs 154 Apr 28 15:24 CBK33ACXX_RNAD8-KORAI1ex4_16s001174-1-1_Serandour_lane7RNAES129D8c14_sequence.txt.gz -> /g/noh/NGS_Library/2016-03-10-CBK33ACXX/lane7/LIB21144_RBA20527/fastq/CBK33ACXX_RNAD8-KORAI1ex4_16s001174-1-1_Serandour_lane7RNAES129D8c14_sequence.txt.gz
lrwxrwxrwx 1 galaxy gbcs 154 Apr 28 15:24 CBK33ACXX_RNAD8-KORAI1ex4_16s001174-1-1_Serandour_lane7RNAES129D8c15_sequence.txt.gz -> /g/noh/NGS_Library/2016-03-10-CBK33ACXX/lane7/LIB21145_RBA20528/fastq/CBK33ACXX_RNAD8-KORAI1ex4_16s001174-1-1_Serandour_lane7RNAES129D8c15_sequence.txt.gz
lrwxrwxrwx 1 galaxy gbcs 154 Apr 28 15:24 CBK33ACXX_RNAD8-KORAI1ex4_16s001174-1-1_Serandour_lane7RNAES129D8c16_sequence.txt.gz -> /g/noh/NGS_Library/2016-03-10-CBK33ACXX/lane7/LIB21140_RBA20523/fastq/CBK33ACXX_RNAD8-KORAI1ex4_16s001174-1-1_Serandour_lane7RNAES129D8c16_sequence.txt.gz
lrwxrwxrwx 1 galaxy gbcs 154 Apr 28 15:24 CBK33ACXX_RNAD8-KORAI1ex4_16s001174-1-1_Serandour_lane7RNAES129D8c18_sequence.txt.gz -> /g/noh/NGS_Library/2016-03-10-CBK33ACXX/lane7/LIB21142_RBA20525/fastq/CBK33ACXX_RNAD8-KORAI1ex4_16s001174-1-1_Serandour_lane7RNAES129D8c18_sequence.txt.gz
(/g/funcgen/galaxy/conda-envs/stocks-client) galaxy@spinoza:/g/funcgen/embase_submissions/kaesper_april2021/rai1/fastq$
```

The STOCKS REST API enables interoperability

- STOCKS REST API available
 - allow integration with other tools e.g. Galaxy and LOOL
 - automate data registration e.g. GC Bridge
 - STOCKS GUI communicates with the server only using the API
 - Write your own client !
- CLI available
 - also uses the REST API
 - extract study information into MAGE-TAB
- Integration with the ITS DMA (soon):
 - Datasets created in STOCKS are pushed to DMA
 - Will also allow import from DMA

Integration with DMA

-  Under dev
-  Existing Routes
-  Potential Routes



Thank you for your attention



Jelle Scholtalbers



Matthias Monfort



Nayeem Reza

Questions ?